

Rat Cd52/CAMPATH-1 antigen ELISA Kit

Catalog No: #EK19483



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cd52/CAMPATH-1 antigen ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat CAMPATH-1 antigen
Target Name	CAMPATH-1 antigen
Other Names	Lymphocyte differentiation antigen B7, CD52, Cdw52, Rb7
Accession No.	Q63064

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.112ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cd52

Note: This product is for in vitro research use only and is not intended for use in humans or animals.