

Rat Aldh5a1/Succinate-semialdehyde dehydrogenase, mitochondrial ELISA Kit

Catalog No: #EK19485

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Aldh5a1/Succinate-semialdehyde dehydrogenase, mitochondrial ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Succinate-semialdehyde dehydrogenase, mitochondrial
Target Name	Succinate-semialdehyde dehydrogenase, mitochondrial
Other Names	Aldehyde dehydrogenase family 5 member A1,NAD(+)-dependent succinic semialdehyde dehydrogenase, Ssadh
Accession No.	P51650

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.056ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Aldh5a1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.