

# Rat Aldh16a1/Aldehyde dehydrogenase family 16 member A1 ELISA Kit

Catalog No: #EK19486



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Aldh16a1/Aldehyde dehydrogenase family 16 member A1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Aldehyde dehydrogenase family 16 member A1
Target Name	Aldehyde dehydrogenase family 16 member A1
Accession No.	Q3T1L0

## Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/ml Sensitivity: 15.68pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Aldh16a1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.