

Mouse Chk1/Serine/threonine-protein kinase Chk1 ELISA Kit

Catalog No: #EK19488



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Chk1/Serine/threonine-protein kinase Chk1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Serine/threonine-protein kinase Chk1
Target Name	Serine/threonine-protein kinase Chk1
Other Names	CHK1 checkpoint homolog,Checkpoint kinase-1,Chk1
Accession No.	O35280

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:43pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Chk1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.