

Mouse Cd69/Early activation antigen CD69 ELISA Kit

Catalog No: #EK19492



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cd69/Early activation antigen CD69 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Early activation antigen CD69
Target Name	Early activation antigen CD69
Other Names	CD69
Accession No.	P37217

Application Details

Assay Type: Sandwich

Detection Range: 78-5000 pg/mL

Sensitivity: less than 39.4 pg/mL

Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul

Assay Time: 3-5h Detection wavelength: 450 nm

gene_full_name: Cd69

Note: This product is for in vitro research use only and is not intended for use in humans or animals.