

Mouse B3galt1/Beta-1,3-galactosyltransferase 1 ELISA Kit

Catalog No: #EK19493



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse B3galt1/Beta-1,3-galactosyltransferase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Beta-1,3-galactosyltransferase 1
Target Name	Beta-1,3-galactosyltransferase 1
Other Names	Beta-1,3-GalTase 1,UDP-Gal:betaGlcNAc beta 1,3-galactosyltransferase-I,UDP-galactose:beta-N-acetyl-glucosamine-beta-1,3-galactosyltransferase 1
Accession No.	O54904

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.162ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:B3galt1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.