

Human ANTXR2/Anthrax toxin receptor 2 ELISA Kit

Catalog No: #EK19494



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ANTXR2/Anthrax toxin receptor 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Anthrax toxin receptor 2
Target Name	Anthrax toxin receptor 2
Other Names	Capillary morphogenesis gene 2 protein,CMG-2,CMG2
Accession No.	P58335

Application Details

Assay Type: Sandwich Detect Range:1.56-100ng/mL Sensitivity:0.78ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:ANTXR2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.