

Mouse Bai3/Brain-specific angiogenesis inhibitor 3 ELISA Kit

Catalog No: #EK19495



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Bai3/Brain-specific angiogenesis inhibitor 3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Brain-specific angiogenesis inhibitor 3
Target Name	Brain-specific angiogenesis inhibitor 3
Accession No.	Q80ZF8

Application Details

Assay Type: Sandwich Detect Range: 78-5000 pg/mL Sensitivity: 39 pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bai3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.