

# Mouse Crhr1/Corticotropin-releasing factor receptor 1 ELISA Kit

Catalog No: #EK19496



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Crhr1/Corticotropin-releasing factor receptor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Corticotropin-releasing factor receptor 1
Target Name	Corticotropin-releasing factor receptor 1
Other Names	CRF-R-1,Corticotropin-releasing hormone receptor 1,CRH-R-1,Crhr
Accession No.	P35347

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.083ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Crhr1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.