

Bovine CLDN3/Claudin-3 ELISA Kit

Catalog No: #EK19498



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine CLDN3/Claudin-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Claudin-3
Target Name	Claudin-3
Accession No.	Q765N9

Application Details

Assay Type: Sandwich Detect Range: 0.625-40ng/mL Sensitivity: 0.38ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CLDN3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.