

# Mouse Cldn3/Claudin-3 ELISA Kit

Catalog No: #EK19499



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cldn3/Claudin-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Claudin-3
Target Name	Claudin-3
Other Names	Clostridium perfringens enterotoxin receptor 2,CPE-R 2,Cpetr2
Accession No.	Q9Z0G9

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.085ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Cldn3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.