

Rat Cldn3/Claudin-3 ELISA Kit

Catalog No: #EK19500



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cldn3/Claudin-3 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Claudin-3
Target Name	Claudin-3
Other Names	Rat ventral prostate.1 protein,RVP1
Accession No.	Q63400

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.12ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cldn3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.