

# Rat Cldn16/Claudin-16 ELISA Kit

Catalog No: #EK19501



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Cldn16/Claudin-16 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Claudin-16
Target Name	Claudin-16
Other Names	Paracellin-1
Accession No.	Q91Y55

## Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.426ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Cldn16

Note: This product is for in vitro research use only and is not intended for use in humans or animals.