

Mouse Clc/Eosinophil lysophospholipase ELISA Kit

Catalog No: #EK19502



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Clc/Eosinophil lysophospholipase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Eosinophil lysophospholipase
Target Name	Eosinophil lysophospholipase
Other Names	Charcot-Leyden crystal protein homolog,CLC,Galectin-10,Gal-10,Lysolecithin acylhydrolase,Clcp
Accession No.	P97400

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Clc

Note: This product is for in vitro research use only and is not intended for use in humans or animals.