

# Mouse Clu/Clusterin ELISA Kit

Catalog No: #EK19508



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Clu/Clusterin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Clusterin
Target Name	Clusterin
Other Names	Apolipoprotein J,Apo-J,Clustrin,Sulfated glycoprotein 2,SGP-2,Apoj,Msgp-2
Accession No.	Q06890

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.063ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Clu

Note: This product is for in vitro research use only and is not intended for use in humans or animals.