

Rat Cndp2/Cytosolic non-specific dipeptidase ELISA Kit

Catalog No: #EK19514



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Rat Cndp2/Cytosolic non-specific dipeptidase ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Rt |
| Specificity | Natural and recombinant Rat Cytosolic non-specific dipeptidase |
| Target Name | Cytosolic non-specific dipeptidase |
| Other Names | CNDP dipeptidase 2 |
| Accession No. | Q6Q0N1 |

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.061ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cndp2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.