

Rat Cndp1/Beta-Ala-His dipeptidase ELISA Kit

Catalog No: #EK19516



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cndp1/Beta-Ala-His dipeptidase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Beta-Ala-His dipeptidase
Target Name	Beta-Ala-His dipeptidase
Other Names	CNDP dipeptidase 1,Carnosine dipeptidase 1,Cn1
Accession No.	Q66HG3

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.081ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cndp1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.