

Rat Asgr1/Asialoglycoprotein receptor 1 ELISA Kit

Catalog No: #EK19520



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Asgr1/Asialoglycoprotein receptor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Asialoglycoprotein receptor 1
Target Name	Asialoglycoprotein receptor 1
Other Names	ASGP-R 1,Hepatic lectin 1,HL-1,Asgr-1
Accession No.	P02706

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.087ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Asgr1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.