

Rat Csk/Tyrosine-protein kinase CSK ELISA Kit

Catalog No: #EK19523



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Csk/Tyrosine-protein kinase CSK ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Tyrosine-protein kinase CSK
Target Name	Tyrosine-protein kinase CSK
Other Names	C-Src kinase
Accession No.	P32577

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.083ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Csk

Note: This product is for in vitro research use only and is not intended for use in humans or animals.