

Mouse Cast/Calpastatin ELISA Kit

Catalog No: #EK19525



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cast/Calpastatin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Calpastatin
Target Name	Calpastatin
Other Names	Calpain inhibitor
Accession No.	P51125

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 26pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cast

Note: This product is for in vitro research use only and is not intended for use in humans or animals.