

Bovine CHGA/Chromogranin-A ELISA Kit

Catalog No: #EK19527



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine CHGA/Chromogranin-A ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Chromogranin-A
Target Name	Chromogranin-A
Other Names	CgA, Pituitary secretory protein I, SP-I
Accession No.	P05059

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.78pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CHGA

Note: This product is for in vitro research use only and is not intended for use in humans or animals.