

Mouse Adora1/Adenosine receptor A1 ELISA Kit

Catalog No: #EK19531



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Adora1/Adenosine receptor A1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Adenosine receptor A1
Target Name	Adenosine receptor A1
Accession No.	Q60612

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.059ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Adora1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.