

Dog CTLA4/Cytotoxic T-lymphocyte protein 4 ELISA Kit

Catalog No: #EK19532



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Dog CTLA4/Cytotoxic T-lymphocyte protein 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Dog
Specificity	Natural and recombinant Dog Cytotoxic T-lymphocyte protein 4
Target Name	Cytotoxic T-lymphocyte protein 4
Other Names	Cytotoxic T-lymphocyte-associated antigen 4,CTLA-4,CD152
Accession No.	Q9XS11

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:39.2pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CTLA4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.