

Bovine CDK4/Cyclin-dependent kinase 4 ELISA Kit

Catalog No: #EK19534



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bovine CDK4/Cyclin-dependent kinase 4 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cyclin-dependent kinase 4
Target Name	Cyclin-dependent kinase 4
Other Names	Cell division protein kinase 4
Accession No.	Q32KY4

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.058ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CDK4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.