

Mouse Cdk4/Cyclin-dependent kinase 4 ELISA Kit

Catalog No: #EK19535



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Cdk4/Cyclin-dependent kinase 4 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Cyclin-dependent kinase 4 |
| Target Name | Cyclin-dependent kinase 4 |
| Other Names | CRK3, Cell division protein kinase 4, PSK-J3, Crk3 |
| Accession No. | P30285 |

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.3pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cdk4

Note: This product is for in vitro research use only and is not intended for use in humans or animals.