

Bovine BPI/Bactericidal permeability-increasing protein ELISA Kit

Catalog No: #EK19537



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Bovine BPI/Bactericidal permeability-increasing protein ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Bovine |
| Specificity | Natural and recombinant Bovine Bactericidal permeability-increasing protein |
| Target Name | Bactericidal permeability-increasing protein |
| Other Names | BPI |
| Accession No. | P17453 |

Application Details

Assay Type: Sandwich Detect Range: 0.78-50ng/mL Sensitivity: 0.432ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: BPI

Note: This product is for in vitro research use only and is not intended for use in humans or animals.