

Mouse Bpi/Bactericidal permeability-increasing protein ELISA Kit

Catalog No: #EK19538



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Bpi/Bactericidal permeability-increasing protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Bactericidal permeability-increasing protein
Target Name	Bactericidal permeability-increasing protein
Other Names	BPI
Accession No.	Q67E05

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.82pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bpi

Note: This product is for in vitro research use only and is not intended for use in humans or animals.