

Rat Cotl1/Coactosin-like protein ELISA Kit

Catalog No: #EK19542



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Cotl1/Coactosin-like protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Coactosin-like protein
Target Name	Coactosin-like protein
Other Names	Clp
Accession No.	B0BNA5

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.089ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cotl1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.