

# Human ASPH/Aspartyl/asparaginyl beta-hydroxylase ELISA Kit

Catalog No: #EK19543



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ASPH/Aspartyl/asparaginyl beta-hydroxylase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Aspartyl/asparaginyl beta-hydroxylase
Target Name	Aspartyl/asparaginyl beta-hydroxylase
Other Names	Aspartate beta-hydroxylase,ASP beta-hydroxylase,Peptide-aspartate beta-dioxygenase,BAH
Accession No.	Q12797

## Application Details

Assay Type: Sandwich Detect Range:78-5000pg/mL Sensitivity:32pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids. Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:ASPH

Note: This product is for in vitro research use only and is not intended for use in humans or animals.