

Mouse Ccl25/C-C motif chemokine 25 ELISA Kit

Catalog No: #EK19547



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Ccl25/C-C motif chemokine 25 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse C-C motif chemokine 25
Target Name	C-C motif chemokine 25
Other Names	Chemokine TECK, Small-inducible cytokine A25, Thymus-expressed chemokine, Scya25, Teck
Accession No.	O35903

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ccl25

Note: This product is for in vitro research use only and is not intended for use in humans or animals.