

Pig CCL25/C-C motif chemokine 25 ELISA Kit

Catalog No: #EK19548



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Pig CCL25/C-C motif chemokine 25 ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Pig |
| Specificity | Natural and recombinant Pig C-C motif chemokine 25 |
| Target Name | C-C motif chemokine 25 |
| Other Names | Chemokine ligand 25, Small-inducible cytokine A25 |
| Accession No. | Q4PR21 |

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CCL25

Note: This product is for in vitro research use only and is not intended for use in humans or animals.