

# Mouse Aif1/Allograft inflammatory factor 1 ELISA Kit

Catalog No: #EK19550



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Aif1/Allograft inflammatory factor 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Allograft inflammatory factor 1
Target Name	Allograft inflammatory factor 1
Other Names	AIF-1, Ionized calcium-binding adapter molecule 1, Iba1
Accession No.	O70200

## Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.053ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: Aif1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.