

Mouse Ar/Androgen receptor ELISA Kit

Catalog No: #EK19552



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Ar/Androgen receptor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Androgen receptor
Target Name	Androgen receptor
Other Names	Dihydrotestosterone receptor, Nuclear receptor subfamily 3 group C member 4, Nr3c4
Accession No.	P19091

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 12pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Ar

Note: This product is for in vitro research use only and is not intended for use in humans or animals.