

Rat Atp2b1/Plasma membrane calcium-transporting ATPase 1 ELISA Kit

Catalog No: #EK19555



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Atp2b1/Plasma membrane calcium-transporting ATPase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Plasma membrane calcium-transporting ATPase 1
Target Name	Plasma membrane calcium-transporting ATPase 1
Other Names	PMCA1, Plasma membrane calcium ATPase isoform 1, Plasma membrane calcium pump isoform 1
Accession No.	P11505

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 28pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids. Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Atp2b1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.