

Mouse Atp2b2/Plasma membrane calcium-transporting ATPase 2 ELISA Kit

Catalog No: #EK19557



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Atp2b2/Plasma membrane calcium-transporting ATPase 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Plasma membrane calcium-transporting ATPase 2
Target Name	Plasma membrane calcium-transporting ATPase 2
Other Names	PMCA2, Plasma membrane calcium ATPase isoform 2, Plasma membrane calcium pump isoform 2, Pmca2
Accession No.	Q9R0K7

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.085ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Atp2b2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.