

Rat Atp2b2/Plasma membrane calcium-transporting ATPase 2 ELISA Kit

Catalog No: #EK19558



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Atp2b2/Plasma membrane calcium-transporting ATPase 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Plasma membrane calcium-transporting ATPase 2
Target Name	Plasma membrane calcium-transporting ATPase 2
Other Names	PMCA2,Plasma membrane calcium ATPase isoform 2,Plasma membrane calcium pump isoform 2
Accession No.	P11506

Application Details

Assay Type: Sandwich Detect Range:0.132-20ng/mL Sensitivity:0.061ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Atp2b2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.