

# Human CDX2/Homeobox protein CDX-2 ELISA Kit

Catalog No: #EK19561



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CDX2/Homeobox protein CDX-2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Homeobox protein CDX-2
Target Name	Homeobox protein CDX-2
Other Names	CDX-3,Caudal-type homeobox protein 2,CDX3
Accession No.	Q99626

## Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.089ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:CDX2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.