

# Rat Ahcy/Adenosylhomocysteinase ELISA Kit

Catalog No: #EK19564



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Ahcy/Adenosylhomocysteinase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Adenosylhomocysteinase
Target Name	Adenosylhomocysteinase
Other Names	AdoHcyase,S-adenosyl-L-homocysteine hydrolase
Accession No.	P10760

## Application Details

Assay Type: Sandwich Detect Range:31.2-2000pg/mL Sensitivity:15.68pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Ahcy

Note: This product is for in vitro research use only and is not intended for use in humans or animals.