

Mouse Btnl2/Butyrophilin-like protein 2 ELISA Kit

Catalog No: #EK19565



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Btnl2/Butyrophilin-like protein 2 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Butyrophilin-like protein 2
Target Name	Butyrophilin-like protein 2
Other Names	Gm315,Ng9
Accession No.	O70355

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.08ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Btnl2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.