

# Mouse Cd46/Membrane cofactor protein ELISA Kit

Catalog No: #EK19566



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Cd46/Membrane cofactor protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Membrane cofactor protein
Target Name	Membrane cofactor protein
Other Names	CD46,Mcp
Accession No.	O88174

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/ml Sensitivity:0.092ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name: Cd46

Note: This product is for in vitro research use only and is not intended for use in humans or animals.