

Mouse Cirbp/Cold-inducible RNA-binding protein ELISA Kit

Catalog No: #EK19568



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Cirbp/Cold-inducible RNA-binding protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cold-inducible RNA-binding protein
Target Name	Cold-inducible RNA-binding protein
Other Names	A18 hnRNP, Glycine-rich RNA-binding protein CIRP, Cirp
Accession No.	P60824

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/mL Sensitivity: 0.086ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cirbp

Note: This product is for in vitro research use only and is not intended for use in humans or animals.