

Mouse Agrp/Agouti-related protein ELISA Kit

Catalog No: #EK19573



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Agrp/Agouti-related protein ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Agouti-related protein
Target Name	Agouti-related protein
Other Names	Agrt,Art
Accession No.	P56473

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.08ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name: Agrp

Note: This product is for in vitro research use only and is not intended for use in humans or animals.