

Human CREB1/Cyclic AMP-responsive element-binding protein 1 ELISA Kit

Catalog No: #EK19578



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CREB1/Cyclic AMP-responsive element-binding protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cyclic AMP-responsive element-binding protein 1
Target Name	Cyclic AMP-responsive element-binding protein 1
Other Names	CREB-1
Accession No.	P16220

Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/mL Sensitivity: 39.8pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CREB1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.