

Rat Acta2/Actin, aortic smooth muscle ELISA Kit

Catalog No: #EK19583



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Rat Acta2/Actin, aortic smooth muscle ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Actin, aortic smooth muscle
Target Name	Actin, aortic smooth muscle
Other Names	Alpha-actin-2,Actsa,Actvs
Accession No.	P62738

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.172ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Acta2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.