

# Bovine BHMT/Betaine--homocysteine S-methyltransferase 1 ELISA Kit

Catalog No: #EK19585



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Bovine BHMT/Betaine--homocysteine S-methyltransferase 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Betaine--homocysteine S-methyltransferase 1
Target Name	Betaine--homocysteine S-methyltransferase 1
Accession No.	Q5I597

## Application Details

Assay Type: Sandwich  
Detect Range: 15.6-1000pg/mL  
Sensitivity: 7.86pg/mL  
Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul  
Assay Time: 3-5h  
Detection wavelength: 450 nm  
gene\_full\_name: BHMT

Note: This product is for in vitro research use only and is not intended for use in humans or animals.