

# Mouse Atf1/Cyclic AMP-dependent transcription factor ATF-1 ELISA Kit

Catalog No: #EK19588



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse Atf1/Cyclic AMP-dependent transcription factor ATF-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cyclic AMP-dependent transcription factor ATF-1
Target Name	Cyclic AMP-dependent transcription factor ATF-1
Other Names	cAMP-dependent transcription factor ATF-1,Activating transcription factor 1,TCR-ATF1
Accession No.	P81269

## Application Details

Assay Type: Sandwich Detect Range:15.6-1000pg/mL Sensitivity:5.2pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Atf1

Note: This product is for in vitro research use only and is not intended for use in humans or animals.