

Mouse Adcy10/Adenylate cyclase type 10 ELISA Kit

Catalog No: #EK19594



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Adcy10/Adenylate cyclase type 10 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Adenylate cyclase type 10
Target Name	Adenylate cyclase type 10
Other Names	Germ cell soluble adenylyl cyclase,sAC,Testicular soluble adenylyl cyclase,Sac,Sacy
Accession No.	Q8C0T9

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.15ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Adcy10

Note: This product is for in vitro research use only and is not intended for use in humans or animals.