

# Rat Alox5/Arachidonate 5-lipoxygenase ELISA Kit

Catalog No: #EK19599



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Alox5/Arachidonate 5-lipoxygenase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Arachidonate 5-lipoxygenase
Target Name	Arachidonate 5-lipoxygenase
Other Names	5-LO
Accession No.	P12527

## Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.194ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene\_full\_name:Alox5

Note: This product is for in vitro research use only and is not intended for use in humans or animals.