

# Human ALOX15B/Arachidonate 15-lipoxygenase B ELISA Kit

Catalog No: #EK19600



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human ALOX15B/Arachidonate 15-lipoxygenase B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Arachidonate 15-lipoxygenase B
Target Name	Arachidonate 15-lipoxygenase B
Other Names	15-LOX-B, 15-lipoxygenase 2, 15-LOX-2, Arachidonate 15-lipoxygenase type II, Linoleate 13-lipoxygenase 15-LOb
Accession No.	O15296

## Application Details

Assay Type: Sandwich Detect Range: 6.25-400 nmol/L Sensitivity: 2.9 nmol/L Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: ALOX15B

Note: This product is for in vitro research use only and is not intended for use in humans or animals.