

Mouse Alox15b/Arachidonate 15-lipoxygenase B ELISA Kit

Catalog No: #EK19601



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Alox15b/Arachidonate 15-lipoxygenase B ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Arachidonate 15-lipoxygenase B
Target Name	Arachidonate 15-lipoxygenase B
Other Names	15-LOX-B, 15-lipoxygenase 2, 15-LOX-2, 8S-lipoxygenase, 8S-LOX, Arachidonate 15-lipoxygenase type II, Alox8
Accession No.	O35936

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Alox15b

Note: This product is for in vitro research use only and is not intended for use in humans or animals.